

PicoGreen Assay

*PicoGreen binds to DNA so handle accordingly: wear appropriate PPE.

*Protect PicoGreen stocks and solutions from light.

Step 1: Calculate # of samples and prepare 1X TE:

_____ of samples + 6 Standards = _____

Amount of 1xTE to prepare:

110 μ L (for picogreen solution) x _____ + 140 μ L (for standard dilutions) = _____

Use 20xTE stock (provided) to prepare appropriate amount of 1xTE from above calculation (or use your own, 1X TE= 10 mM Tris 8, 1 mM EDTA).

Remove 140 μ L for use of making Standards in next step.

Step 2: Prepare Standards:

Use 2 μ L of stock λ -DNA of 100 ng/ μ L (provided) + 18 μ L of 1X TE = 20 μ L of 10 ng/ μ L working solution, this is standard #1. Make five serial dilutions using 2 ul working solution into 18 ul TE to make the following standards:

$$\#2) \frac{(2 \mu\text{L})(10\text{ng}/\mu\text{L})}{(20\mu\text{L})} = 1 \text{ ng}/\mu\text{L}$$

$$\#3) \frac{(2\mu\text{L})(1\text{ng}/\mu\text{L})}{(20\mu\text{L})} = 0.1 \text{ ng}/\mu\text{L}$$

$$\#4) \frac{(2\mu\text{L})(0.1\text{ng}/\mu\text{L})}{(20\mu\text{L})} = 0.01 \text{ ng}/\mu\text{L}$$

$$\#5) \frac{(2\mu\text{L})(0.01\text{ng}/\mu\text{L})}{(20\mu\text{L})} = 0.001 \text{ ng}/\mu\text{L}$$

Put 10 μ L of the original 10 ng/ μ L λ -DNA standard into a designated well on your plate (typically a 96 well full skirted black plate from MJ Research), and continue by adding 10 ul of dilutions 2-5 into the next series of wells. Put 10 ul TE with no DNA into the next well. Now you should have 6 standard wells containing 10, 1, 0.1, 0.01, 0.001, and 0 ng/ μ L of λ -DNA (make note of where each was put).

Step 3: Prepare PicoGreen Solution:

Dilute the 200X PicoGreen stock to 1X PicoGreen with 1X TE prepared above (excluding amount removed for Standards step).

Put 100 μ L of 1xPicoGreen in each well containing 10 μ L of Standard with pipetting up and down to mix.

Continue to add 100 μ L of 1xPicoGreen to each well that will be used to measure your samples.

Add 10 μ L of each sample to wells that have 100 μ L of 1X PicoGreen and pump up and down to mix.

Notes: Your standard curve should be linear to the 0.01 ng/ μ L sample, then it will bottom out, and 0.001 ng/ μ L will probably be the same as 0. Readings can be done right away, but they don't have to be: the dye:DNA mixes are stable for at least 24 hrs at 4°. It may also be possible to do readings in a 384 well format, talk to us if you're interested in trying that out.