




**Amplicon Sequencing**




emPCR allows for single molecule sequencing without cloning the molecules of the PCR sample sequences into bacteria. The emPCR and sequencing reactions are performed following minor variations on the standard Genome Sequencer System procedures.




**Amplicon Sequencing**




For data analysis, a complete software solution is available. The software provided is specifically tailored for the determination and quantification of sequence deviations. It makes available very powerful tools and a convenient interface able to identify rare variants, quantify their prevalence, and even verify whether variations observed at different sequencing positions are linked to each other, by looking at individual reads.



**Amplicon Sequencing**



During the experimental design phase, **fusion primers** are designed for the generation of PCR products. **These PCR products, resulting from a simple PCR amplification step,** represent the Amplicon DNA library that contains the genomic spots that interest you. At the end of this experimental step, you will have an **Amplicon DNA library ready** for the amplification phase of the Genome Sequencer System.



**Fusion Primer**

Fixed sequencing primer: 19-mer binds to Capture Beads

Target-specific sequencing primer: anneals amplification & sequencing primers

includes sequencing key

5' GCCTCCCTCGCGCATCAG 3'

5' GCCTTGCCAGCCCGTCAG 3'

The primers used to generate amplicon libraries are composed of two sequence parts: a **fixed sequencing primer** and a **target-specific sequencing primer**. The fixed sequencing primer at the 5' end is a Genome Sequencer System-specific **19-mer**. This primer **binds to the Capture Beads, anneals the emPCR amplification primers** and the sequencing primer, and ends with the **sequencing key "TCAG"**.

**Amplicon Sequencing**

Experimental Design → emPCR → Sequencing Run

**Fusion Primer**

Fixed sequencing primer | Fusion Primer | Target-specific sequencing primer

Primer A  
5' GCCTCCCTCGCGCCATCAG 3'

Primer B  
5' GCCTTGCCAGCCCGCTCAG 3'

GS emPCR Kit III (Amplicon B)

Primer A | Primer B

Primer B is complementary to the capture oligonucleotide which is covalently linked to the DNA Capture Beads in the GS emPCR Kit II. This kit allows you to sequence the amplified fragment from Primer A.

**Amplicon Sequencing**

Experimental Design → emPCR → Sequencing Run

**Fusion Primer**

Fixed sequencing primer | Fusion Primer | Target-specific sequencing primer

> user designed, target-specific sequence  
> 20-25 bp

target


The **target-specific sequencing primer at the 3' end** is typically **20-25 base pairs** in length and designed to anneal to **either side** of the target to be sequenced. This requires some **knowledge of the target sequence**, and works best on target sequences that are not too variable.

**Amplicon Sequencing**

Experimental Design → Inland cloning ampPCR → Sequencing Run

**Amplicon**

PCR Amplification




The current GS technology generates reads at an average length of greater than 200 base pairs. The distance between the two sequence-specific primers must be carefully considered if the "inserts" are to be sequenced in their entirety.

In general **we recommend bidirectional sequencing** to achieve high accuracy.

PCR Amplification

> Unidirectional sequencing



Amplicon length  $\leq 400$  bp

If unidirectional sequencing without any overlap is sufficient, the amplicon should be no longer than about 400 base pairs from key to key


**Amplicon Sequencing**

Experimental Design → Inland cloning ampPCR → Sequencing Run

**Amplicon**

PCR Amplification


> Unidirectional sequencing



For example, if the target-specific components of the Fusion Primers are each 21-mers, you will read about 180 nucleotides into the insert from each side - thus covering a 360 basepair insert.



PCR Amplification



> Bidirectional sequencing



Amplicon length: 200 bp

If bidirectional sequencing is required, the amplicon should be **no longer than about 200 basepairs from key to key**. With the same example of **21 basepair target-specific components of the Fusion Primers**, you can read a **160 basepair insert from both ends**, so the amplicon can be this length plus the **target-specific primer** and the **key** on either side.







Amplicon libraries can be read from either end, or from both the "A" and "B" primer ends.

To read from Primer A, process the Amplicon library using the GS emPCR Kit II; to read from Primer B, process the library using the GS emPCR Kit III.

To read from both ends, **process an aliquot of the library using each of the GS emPCR Kits II and III and pool the two clonally amplified libraries before loading them onto the PicoTiterPlate device** for sequencing.

Amplicon Sequencing

Experimental Design → Instand cloning emPCR → Sequencing Run → Data Analysis

Sequencing-by-synthesis

PicoTiterPlate device

Sequencing-by-synthesis

(Bead depositing)

(Bead depositing)


(Sequencing process)

→ Generation of a chemiluminescent signal

Another special feature of the Genome Sequencer System is its sequencing technology. Sequencing-by-synthesis is performed **simultaneously** in the **open wells** of a fiber-optic slide called a **PicoTiterPlate device**. For Amplicon Sequencing, and for all the applications of the Genome Sequencer System, sequencing-by-synthesis follows the same procedures.

First, the **DNA beads are loaded** together with **enzyme** and **packing beads** into the wells of the PicoTiterPlate device.

When the sequencing process starts, nucleotides are delivered to the wells in a **cyclically** repeating flow sequence. The sequencing-by-synthesis method is based on the generation of a chemiluminescent signal. The signal intensity at each nucleotide flow, for a particular well, indicates the number of nucleotides - if any - that were incorporated. If no nucleotide is incorporated during a nucleotide flow, no light signal is generated, and the complementary strand is not elongated.



**Amplicon Variant Analyzer (AVA) Software**

Overview | Project | Computations | Variants | Global Align | Consensus Align | **Flowgrams**

**Project**

Name: EGFR

Location: /home/training/ampProjects/trainingProject001

Description: Study of EGPR somatic mutations

Survey of exons 18-22

This pilot study will guide the selection of future samples and facilitate the discovery of

**Summary**

References: 5

Amplification: 11

Read Data: 7

Variants: 4



## GS FLX Titanium



## GS Junior System



<b>Throughput</b>	400-600 million high-quality, filter-passed bases per run* 1 billion bases per day
<b>Run Time</b>	10 hours
<b>Read Length</b>	Average length = 400 bases
<b>Accuracy</b>	Q20 read length of 400 bases (99% at 400 bases and higher for prior bases)
<b>Reads per run</b>	>1 million high-quality reads
<b>Data</b>	Trace data accepted by NCBI since 2005
<b>Computing Requirements</b>	Cluster recommended (Roche GS FLX Titanium Cluster available)
<b>Robustness</b>	No complex optics or lasers; reagents have long shelf life

### System Performance

<b>Throughput</b>	35 million high-quality, filtered bases per run*
<b>Run Time</b>	10 hours sequencing 2 hours data processing
<b>Avg. Read Length</b>	400 bases*
<b>Accuracy</b>	Q20 read length of 400 bases (99% accuracy at 400 bases)
<b>Reads per Run</b>	100,000 shotgun; 70,000 amplicon
<b>Sample Input</b>	gDNA, amplicons, cDNA, or BACs depending on the application
<b>Physical Dimensions</b>	40 cm wide x 60 cm deep x 40 cm high (the size of a laser printer) Weight = 55 lbs.
<b>Computing</b>	Linux-based OS on HP desktop computer included. All software is point-and-click.

\*Typical results. Average read length and number of reads depend on specific sample and genomic characteristics